

Technical Data Sheet

Eye Wash Station (EPS) complete



Order no.: 277 999

EAN: 4013207277997

Product designation:

Eye Wash Station (EPS) complete, with 2 EKASTU Eye Wash Bottles

consisting of EKASTU Eye Wash Bottle ECO with funnel, approx. 600 ml FD and EKASTU Eye Wash Bottle MINI-ECO with funnel, approx. 200 ml FD and one Eye Wash Station for protected storage.

EKASTU Eye Wash Bottles are manufactured in accordance with DIN EN 15154-4 and are a Class I Medical Device, in accordance with the Regulation (EU) 2017/745.



DESCRIPTION

Eye Wash Station (EPS) complete including EKASTU Eye Wash Bottle ECO with funnel, approx. 600 ml FD and EKASTU Eye Wash Bottle MINI-ECO with funnel, approx. 200 ml FD.

Eye Wash Bottles according to DIN EN 15154-4

Protects Eye Wash Bottles against coldness and warmth and the water content in the Eye Wash Bottle has a temperature delta of approx. 38 °C (corresponds to +/- 19 °C). The eye wash station protects the contents of the Eye Wash Bottles almost 5 times better against warmth and coldness

For immediate availability in hazardous areas, easy-to-open lid, dust-proof and easy wall mounting

MAIN FEATURES



MATERIALS

- Expandable polystyrene (EPS)
- EKASTU Eye Wash Bottle ECO with funnel, approx. 600 ml FD see technical data sheet for Item no. 177 111
- EKASTU Eye Wash Bottle MINI-ECO with funnel, approx. 200 ml FD see technical data sheet for Item no. 177 112

The materials used are particularly kind of the skin and therefore have no negative effect on human health.

WEIGHT

- approx. 1190 g

MARKING

Eye Wash Bottle

- Trade mark of manufacturer (LOGO)
- Manufacturer/- address
- CH-REP
- Product designation
- Applied standard
- Consider the information brochure
- LOT-No.
- Expiry of shelf-life (MM/YYYY)
- CE-Identification
- Storage conditions

DONNING AND USING

- Protect from sunlight, chemicals, heat and cold. Eye Wash Station (EPS) must be visually inspected at regular intervals
- Do not open the Eye Wash Bottle until immediately before use
- Protect from sunlight, chemicals, heat and cold. Observe storage conditions on the packaging. The eyewash bottle must be visually inspected at regular intervals

WARNING REFERENCE

- Use in accordance with the requirements only!
- Please consider the information brochure!

CERTIFICATIONS

- EKASTU Eye Wash Bottles are manufactured in accordance with DIN EN 15154-4 and are a Class I Medical Device, in accordance with the Regulation (EU) 2017/745

USER REFERENCES

EKASTU Safety GmbH guarantees the indicated achievement according to class and type. It is to be noted that the laboratory test values can considerably deviate from those, which are reached in practice. The user must read and understand all functional information. Use the Eye Wash Bottle with this information brochure, the relevant valid statutory regulations and the safety requirements of the profession associations for example DGUV Information 204-022 "First Aid in Operation".

TD 277999 English_01-10.24