



High Tech Engineering

LabHEAT® residual current detection Type KM-PSA16 and KM-PSS16

KM-PSA16

RCD (residual current detection) adapter

Adapter for personal protection, the operation of heating mantles with RCD is imperative!

This adapter is to switch between heating mantle and source of power socket.

Technical details

- tripping current > 30 mA
- tripping time < 30 ms
- max. power load 3600 W (16A)
- system of protection IP40
- nominal voltage 230 V AC
- GS-tested



KM-PSS16

RCD (residual current detection) plug

Technical identical to KM-PSA16, but this RCD is to install fixed on the power supply cable instead the European schuko plug.



Ordering details

Typ		Art.-Nr.
KM-PSA16	RCD-adapter	97150200
KM-PSS16	RCD-plug	97150300

