




Fig. Test chamber Art. 7923506, complete with Cover plate 5192201 and accessories

- 1 7923506 **All-glass test chamber. L 360* B.230*H260**
Edge ground, Clear glass ca. L. 360* B.230* H.260 Contents ca. 20 Liter

For example, used as corrosion - Test Chamber
The chamber is made of seamless glass

Options: Cover plate D. 8 mm 450 *320 mm
Art. 5192201 with added Duran glass fittings: 1* GL 18 (Hole drilling 12 mm)
1* Sleeve NS 29/32 (Hole drilling 30 mm)

If desired, other valves or holes in the lid possible

- 2 5192201 **Cover plate D. 8 mm 390 *260 mm**
with added Duran glass fittings:
1* GL 18 (Hole 15 mm)
1* Sleeve NS 29/32 (Hole 30 mm)
Fittings with a broad base for gluing.

Option:

- 3 7923506P **Fine ground edge**
For high density, extra-flat ground for Test Chamber. Art 7923506P. The Chamber is additionally ground to reach a dense cover edition of the reach.

7923506-Ganzglas-Testkammer-eng