



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	802CE.A
GENERAL DESCRIPTION	MINI URISPONGE IN BULK AUTOMATION 300PCS
DETAILED DESCRIPTION	1 box containing 6 envelopes each containing 50 tubes with yellow screw cap with a stick with strung sponges.

PRODUCT COMPONENTS

SAMPLING TUBE	TUBE	SPECIFICATIONS
		TOTAL LENGTH: 81.30 mm $\pm$ 0.20 mm OUTER TUBE DIAMETER: 13.00 mm $\pm$ 0.20 mm INNER TUBE DIAMETER: 11.50 mm $\pm$ 0.20 mm MATERIAL: POLYPROPYLENE
CAP	CAP YELLOW	SPECIFICATIONS
		CAP HEIGHT: 14.95 mm $\pm$ 0.25 mm LOWER CAP DIAMETER: 16.80 mm $\pm$ 0.20 mm HIGHER CAP DIAMETER: 16.00 mm $\pm$ 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE

This bulletin (Rev.07-26/07/2023) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

SPONGE	SPONGE	SPECIFICATIONS
		SPONGE DIAMETER: 8.0 mm $\pm$ 0.5 mm LENGTH OF SPONGE: 11.0 mm $\pm$ 1.0 mm MATERIAL: POLYURETHANE
APPLICATOR	MINI STICK	SPECIFICATIONS
		TOTAL LENGTH: 60.80 mm $\pm$ 2.00 mm SPONGES AREA: 24.00 mm $\pm$ 1.00 mm DIAMETER AT CAP CONNECTION: 2.50 mm $\pm$ 0.10 mm MATERIAL: POLYSTYRENE

PACKAGING SPECIFICATION

PACKAGING	BULK
BOX QTY	300
ENVELOPE QTY	50

IMPORTANT INFO

STERILIZATION METHOD	IONIZING IRRADIATION
STORAGE TEMPERATURE	2/25°C-DO NOT FREEZE
SHELF LIFE	12 MONTHS

This bulletin (Rev.07-26/07/2023) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.