

PRODUCT DATA SHEET



Product description

Catalogue Number	5-203-C4-0
Product description	AHN myTip® LFT 1000 µL, clear, sterile, racked
Sales Unit	case / 4 x 8 x 96 tips

Product characteristics

Graduations	yes
Filter	yes
Filter Specifications	18-40 µm, PE, hydrophob
Low Retention	yes
Sterilised	yes
Autoclavable (max. °C, period)	no
Purity	Free of RNase, DNase, Human DNA, Pyrogens, PCR Inhibitors

Size

Max. volume	1000 µL
Length of product	84.30 mm

Material/colours

Product material	polypropylene (PP)
Colour of product	clear

Packaging

Type of Packaging	Racked
Pcs./Rackbox	96 tips
Pcs./Pack	768 tips
Pcs./Case	3072 tips
Dimension Case Carton (H x W x D)	399 x 273 x 400 mm

Storage and Handling

Temperature	below 35°C
UV light	protect from UV light

Note: Improper storage can initiate degradation, which results in colour changes and can have negative effects on the physical properties of the products.

Safety: The product is not classified as dangerous.

Recycling: The product is recyclable.

This document was sent online and is valid without signature

CEOs

Alexander Spector
Chirag Shah
Henrique Hosken
Tejas Shah
Torsten Buchwald

Company Address

AHN Biotechnologie GmbH
Uthleber Weg 14
D-99734 Nordhausen

Contact

T: +49 (0) 3631 65242-0
F: +49 (0) 3631 65242-90
E: info@cappahn.com
W: www.cappahn.com

Bank Account

Commerzbank AG €/€
Konto-Nr.: 6065155
BLZ: 82040000
BIC: COBADEFF820
IBAN: DE16820400000606515500

Commercial Register

HRB 404781
Amtsgericht Jena
UST-IdNr. DE812717024

