

Product Data Sheet



Product Description

Catalogue number	5030061C
Product description	ExpellPlus 10/20µL XL, extra narrow, pre-sterile, w/filter, hinged racks
Sales Unit	Case / 5x10x96 tips

Product Characteristics

Graduation	Yes
Filter	Yes
Filter specifications	18-40 µm, PE, hydrophob
Low retention	Yes
Sterilised	Yes
Autoclavable	No
Purity	Free of RNase, DNase, Human DNA, Pyrogens, PCR Inhibitors

Size

Max. volume	20 µL
Length of product	45.80mm

Material / Colours

Product material	Polypropylene (PP)
Colour of product	Clear

Packaging

Type of packaging	Racked
Pieces/rackbox	96 tips
Pieces / bag	
Pieces / pack	960 tips
Pieces / case	4800 tips
Dimensions case / carton	296 x 338 x 400 mm

Storage and Handling

Temperature	Below 35°C
UV light	Protect from UV light

Note: *Improper storage can initiate degradation, which results in colour changes and can have negative effects on the physical properties of the products.*

Safety: *The product is not classified as dangerous*

Recycling: *The product is recyclable.*