

CAPP Aero 96

MULTICHANNEL PIPETTES

Easy and light tip mounting and ejection

Excessive use of the thumb is a major cause of repetitive strain injuries and user discomfort in multichannel pipetting.

The unique design of the **CAPP Aero 96 Multichannel Pipette** considerably reduces the thumb workload. Ejection is simply a matter of applying light pressure to the levered ejector arm.

Thanks to the double O-rings on the aluminium and metal nose cones, the **CAPP Aero 96 Multichannel Pipette** eliminates the need to 'hammer' your pipette into the tip rack to mount your tip correctly. You only need to apply a slight even pressure to perfectly secure the tips to the nose cones.

CAPP Aero 96 Multichannel Pipettes are one of the lightest pipettes on the market and produce a consistent and reliable liquid transfer result in one ergonomic solution.

Single channel accuracy and precision in a multichannel package

The **CAPP Aero 96 Multichannel Pipette** achieves single channel accuracy and precision across all channels due to the special design and construction of the **CAPP** piston assembly.

The robust, fully autoclavable pipette body offers easy maintenance, while the unique **CAPP** design, with individually replaceable tip bases, reduces the cost of what were once inevitable repairs.

The **CAPP Aero 96 Multichannel Pipette's** use of anodised aluminium and metal nose cones instead of plastic ones substantially reduces breakage caused by accidental dropping or inaccurate tip mounting.

The majority of internal construction parts are made of lightweight metal components and, along with stainless steel pistons, ensure a highly reliable performance and stability over time in terms of accuracy and precision. This smart, highly ergonomic pipette is built to last.



Ordering Information

CAPP Aero 96 Multichannel Pipettes

Cat. No.	Description	Inaccuracy %	Imprecision %	Colour
C10-8	CAPPAero 96 multi pipettes, 8-channel, 0.5-10µL	5.00 / 2.00	3.60 / 1.00	○
C20-8	CAPPAero 96 multi pipettes, 8-channel, 2-20µL	5.00 / 2.00	3.60 / 1.00	●
C50-8	CAPPAero 96 multi pipettes, 8-channel, 5-50µL	3.20 / 1.40	1.60 / 0.60	●
C100-8	CAPPAero 96 multi pipettes, 8-channel, 10-100µL	2.50 / 1.30	1.20 / 0.50	●
C200-8	CAPPAero 96 multi pipettes, 8-channel, 20-200µL	1.80 / 1.20	0.80 / 0.40	●
C300-8	CAPPAero 96 multi pipettes, 8-channel, 30-300µL	1.66 / 1.40	0.60 / 0.30	●
C10-12	CAPPAero 96 multi pipettes, 12-channel, 0.5-10µL	5.00 / 2.00	3.60 / 1.00	○
C20-12	CAPPAero 96 multi pipettes, 12-channel, 2-20µL	5.00 / 2.00	3.60 / 1.00	●
C50-12	CAPPAero 96 multi pipettes, 12-channel, 5-50µL	3.20 / 1.40	1.60 / 0.60	●
C100-12	CAPPAero 96 multi pipettes, 12-channel, 10-100µL	2.50 / 1.30	1.20 / 0.50	●
C200-12	CAPPAero 96 multi pipettes, 12-channel, 20-200µL	1.80 / 1.20	0.80 / 0.40	●
C300-12	CAPPAero 96 multi pipettes, 12-channel, 30-300µL	1.66 / 1.40	0.60 / 0.30	●



The new CAPP Sharp Reagent Reservoir.

