

# CAPP Rondo

3 D SHAKER

**CAPP Rondo 3D Shaker** operates at a variable speed range from 10 to 120 RPM. The brushless DC motor ensures stable, quiet and maintenance-free operation, which is perfectly suited for extended runs. The three-dimensional rocking and tumbling motion provide thorough homogenous mixing.

The user can save up to nine individual protocols, or use the pulse mode with *Run-Pause-Run* operation at a programmable time interval ranging from 30 to 99 seconds.

**CAPP Rondo 3D Shaker** can be programmed to run from 1 up to 999 minutes, or continuously in the *Infinite Mode*. To assure total stability on the bench it is equipped with non-slip rubber feet and comes with a spiked rubber mat. Optionally, a sticky pad mat can be ordered.

**CAPP Rondo 3D Shaker** is easy to clean and complies with Protection Class IP 21.



## Specifications:

|                        |                                                         |
|------------------------|---------------------------------------------------------|
| Load Capacity          | 5kg                                                     |
| Speed Range            | 10-120 RPM                                              |
| Motor Type             | BLDC (Brushless DC motor)                               |
| Movement               | Tumbling 7°                                             |
| Run Time               | 1-999 minutes or infinite                               |
| Programmable Mode      | Up to 9 user-defined protocols                          |
| Pulse Mode             | <i>Run-Pause-Run</i> programmable from 30 to 99 secs.   |
| Display                | 7 Segment LED                                           |
| Accessories            | Spiked rubber mat (included), Sticky pad mat (optional) |
| Plate size             | 350 x 350mm                                             |
| Protection Class       | IP 21                                                   |
| Dimensions (W x D x H) | 290 x 370 x 200mm                                       |

## Ordering Information

| Cat. No. | Description                      |
|----------|----------------------------------|
| CR3-08X  | CAPPRondo 3D Shaker 5kg capacity |

