

# Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 14.02.2024

Version number 1

Revision: 14.02.2024

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### · 1.1 Product identifier

- Trade name: **AG Refractometer Calibration Fluid**
- Article number: 90-204

### · 1.2 Relevant identified uses of the substance or mixture and uses advised against:

- Product category: PC21 Laboratory chemicals
- Process category: PROC15 Use as laboratory reagent
- Application of the substance / the preparation: Calibration fluid for refractometers

### · 1.3 Details of the supplier of the safety data sheet

#### · Manufacturer/Supplier:

Xylem Analytics Germany GmbH  
Bellingham + Stanley  
Am Achalaich 11  
82362 Weilheim  
Germany  
Tel. +49 881 183-0

- Further information obtainable from: Email: techsupport.bs.uk@xylem.com

### · 1.4 Emergency telephone number: Chemtrec: (USA & Canada) 800-424-9300 (International) 001 703-527-3887

## SECTION 2: Hazards identification

### · 2.1 Classification of the substance or mixture

- Classification according to Regulation (EC) No 1272/2008: The product is not classified, according to the CLP regulation.

### · 2.2 Label elements:

- Labelling according to Regulation (EC) No 1272/2008: Void
  - Hazard pictograms: Void
  - Signal word: Void
  - Hazard statements: Void

### · 2.3 Other hazards No further relevant information available.

## SECTION 3: Composition/information on ingredients

### · 3.2 Mixtures

#### · Description:

Mixture, consisting of the following components:

- 1,2-Propanediol (CAS 57-55-6): <65%
- Water (CAS 7732-18-5): ≥35%

#### · Dangerous components: Void

- Additional information: For the wording of the listed hazard phrases refer to section 16.

## SECTION 4: First aid measures

### · 4.1 Description of first aid measures

- After inhalation: Supply fresh air.

(Contd. on page 2)

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 14.02.2024

Version number 1

Revision: 14.02.2024

**Trade name: AG Refractometer Calibration Fluid**

(Contd. of page 1)

- **After skin contact:**
  - Rinse skin with water [or shower].
  - Take off contaminated clothing.
  - If skin irritation or rash occurs: Get medical advice/attention.
- **After eye contact:** Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- **After swallowing:**
  - Rinse out mouth and then drink plenty of water.
  - If symptoms persist consult doctor.
- **4.2 Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **4.3 Indication of any immediate medical attention and special treatment needed** No further relevant information available.

## SECTION 5: Firefighting measures

- **5.1 Extinguishing media**
  - **Suitable extinguishing agents:** The product is not flammable. Extinguishing agent to suit environment.
  - **For safety reasons unsuitable extinguishing agents:**
    - There are no extinguishing agent restrictions for this substance/mixture.
- **5.2 Special hazards arising from the substance or mixture** No further relevant information available.
- **5.3 Advice for firefighters**
  - **Protective equipment:** No special measures required.
  - **Additional information** Dispose of fire debris and contaminated fire fighting water in accordance with official regulations.

## SECTION 6: Accidental release measures

- **6.1 Personal precautions, protective equipment and emergency procedures**
  - No special measurements required when handled properly.
- **6.2 Environmental precautions:**
  - Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
- **6.3 Methods and material for containment and cleaning up:**
  - Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
  - Wash off residuals with water.
- **6.4 Reference to other sections**
  - See Section 7 for information on safe handling.
  - See Section 8 for information on personal protection equipment.
  - See Section 13 for disposal information.

## SECTION 7: Handling and storage

- **7.1 Precautions for safe handling** No special precautions are necessary if used correctly.
  - **Information about fire - and explosion protection:** No special measures required.
- **7.2 Conditions for safe storage, including any incompatibilities**
  - **Requirements to be met by storerooms and receptacles:** No special requirements.
  - **Information about storage in one common storage facility:** Not required.
  - **Further information about storage conditions:** Store tightly sealed at temperatures between 15 °C and 25 °C.
- **7.3 Specific end use(s)** No further relevant information available.

EU  
(Contd. on page 3)

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 14.02.2024

Version number 1

Revision: 14.02.2024

Trade name: AG Refractometer Calibration Fluid

(Contd. of page 2)

### SECTION 8: Exposure controls/personal protection

#### · 8.1 Control parameters

##### · Ingredients with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

· **Additional information:** The lists valid during the making were used as basis.

#### · 8.2 Exposure controls

· **Appropriate engineering controls** No further data; see section 7.

##### · General protective and hygienic measures:

Avoid contact with the eyes and skin.

Do not inhale gases / fumes / aerosols.

Take off contaminated clothing and wash it before reuse.

Wash hands before breaks and at the end of work.

##### · Individual protection measures, such as personal protective equipment

· **Respiratory protection:** Not required.

· **Hand protection** Protective gloves

· **Material of gloves**

Nitrile rubber, NBR

Recommended thickness of the material:  $\geq 0.12$  mm

· **Penetration time of glove material** 60-60 min

· **Eye/face protection** Safety glasses

##### · Environmental exposure controls

Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

### SECTION 9: Physical and chemical properties

#### · 9.1 Information on basic physical and chemical properties

##### · General Information

· **Physical state**

Fluid

· **Colour:**

Colourless

· **Odour:**

Odourless

· **Melting point/freezing point:**

Undetermined.

· **Boiling point or initial boiling point and boiling range**

Undetermined.

· **Flammability**

Product is not flammable.

· **Flash point:**

>99 °C (c.c.)

· **Decomposition temperature:**

Not determined.

· **pH (100 g/l) at 20 °C**

6-8

· **Viscosity:**

· **Kinematic viscosity**

Not determined.

· **Dynamic:**

Not determined.

· **Solubility**

· **water:**

Fully miscible.

· **Partition coefficient n-octanol/water (log value)**

Not determined.

· **Vapour pressure:**

Not determined.

· **Density and/or relative density**

· **Density at 20 °C:**

1 g/cm<sup>3</sup>

Not determined.

· **Relative density**

Not determined.

· **Vapour density**

Not determined.

(Contd. on page 4)

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 14.02.2024

Version number 1

Revision: 14.02.2024

Trade name: AG Refractometer Calibration Fluid

(Contd. of page 3)

|                                                                                        |                                               |
|----------------------------------------------------------------------------------------|-----------------------------------------------|
| · <b>9.2 Other information</b>                                                         |                                               |
| · <b>Appearance:</b>                                                                   |                                               |
| · <b>Form:</b>                                                                         | Fluid                                         |
| · <b>Important information on protection of health and environment, and on safety.</b> |                                               |
| · <b>Ignition temperature:</b>                                                         | Product is not selfigniting.                  |
| · <b>Explosive properties:</b>                                                         | Product does not present an explosion hazard. |
| · <b>Change in condition</b>                                                           |                                               |
| · <b>Evaporation rate</b>                                                              | Not determined.                               |
| · <b>Information with regard to physical hazard classes</b>                            |                                               |
| · <b>Explosives</b>                                                                    | Void                                          |
| · <b>Flammable gases</b>                                                               | Void                                          |
| · <b>Aerosols</b>                                                                      | Void                                          |
| · <b>Oxidising gases</b>                                                               | Void                                          |
| · <b>Gases under pressure</b>                                                          | Void                                          |
| · <b>Flammable liquids</b>                                                             | Void                                          |
| · <b>Flammable solids</b>                                                              | Void                                          |
| · <b>Self-reactive substances and mixtures</b>                                         | Void                                          |
| · <b>Pyrophoric liquids</b>                                                            | Void                                          |
| · <b>Pyrophoric solids</b>                                                             | Void                                          |
| · <b>Self-heating substances and mixtures</b>                                          | Void                                          |
| · <b>Substances and mixtures, which emit flammable gases</b>                           |                                               |
| in contact with water                                                                  | Void                                          |
| · <b>Oxidising liquids</b>                                                             | Void                                          |
| · <b>Oxidising solids</b>                                                              | Void                                          |
| · <b>Organic peroxides</b>                                                             | Void                                          |
| · <b>Corrosive to metals</b>                                                           | Void                                          |
| · <b>Desensitised explosives</b>                                                       | Void                                          |

## SECTION 10: Stability and reactivity

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability**
  - **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **10.3 Possibility of hazardous reactions** No dangerous reactions known.
- **10.4 Conditions to avoid** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** No dangerous decomposition products known.

## SECTION 11: Toxicological information

- **11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008**
  - **Acute toxicity**  
No quantitative toxicity data are available for this product.  
Based on available data, the classification criteria are not met.
  - **Skin corrosion/irritation** Based on available data, the classification criteria are not met.
  - **Serious eye damage/irritation** Based on available data, the classification criteria are not met.
  - **Respiratory or skin sensitisation** Based on available data, the classification criteria are not met.
  - **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
  - **Carcinogenicity** Based on available data, the classification criteria are not met.
  - **Reproductive toxicity** Based on available data, the classification criteria are not met.

(Contd. on page 5)

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 14.02.2024

Version number 1

Revision: 14.02.2024

Trade name: AG Refractometer Calibration Fluid

(Contd. of page 4)

- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.

### 11.2 Information on other hazards

#### · Endocrine disrupting properties

None of the ingredients is listed.

## SECTION 12: Ecological information

### 12.1 Toxicity

- **Aquatic toxicity:** No further relevant information available.

### 12.2 Persistence and degradability

No further relevant information available.

### 12.3 Bioaccumulative potential

No further relevant information available.

### 12.4 Mobility in soil

No further relevant information available.

### 12.5 Results of PBT and vPvB assessment

Not applicable.

### 12.6 Endocrine disrupting properties

The product does not contain substances with endocrine disrupting properties.

### 12.7 Other adverse effects

#### · General notes:

Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

#### · Additional ecological information:

- **General notes:** Not known to be hazardous to water.

#### · Additional ecological information:

Not known to be hazardous to water.

## SECTION 13: Disposal considerations

### 13.1 Waste treatment methods

#### · Recommendation

Disposal must comply with the relevant local regulations. We recommend that you contact either the authorities in charge or approved waste disposal companies which will advise you on how to dispose the special waste.

#### · Uncleaned packaging:

##### · Recommendation:

Empty contaminated packagings thoroughly. They may be recycled after thorough and proper cleaning.

Packagings that may not be cleansed are to be disposed of in the same manner as the product.

- **Recommended cleansing agents:** Water, if necessary together with cleansing agents.

## SECTION 14: Transport information

### 14.1 UN number or ID number

- **ADR/RID, IMDG, IATA** Void

### 14.2 UN proper shipping name

- **ADR/RID, IMDG, IATA** Void

### 14.3 Transport hazard class(es)

- **ADR/RID, ADN, IMDG, IATA**
- **Class** Void

### 14.4 Packing group

- **ADR/RID, IMDG, IATA** Void

(Contd. on page 6)

# Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 14.02.2024

Version number 1

Revision: 14.02.2024

**Trade name: AG Refractometer Calibration Fluid**

(Contd. of page 5)

|                                                                       |                                                      |
|-----------------------------------------------------------------------|------------------------------------------------------|
| · <b>14.5 Environmental hazards:</b>                                  | Not applicable.                                      |
| · <b>14.6 Special precautions for user</b>                            | Not applicable.                                      |
| · <b>14.7 Maritime transport in bulk according to IMO instruments</b> | Not applicable.                                      |
| · <b>Transport/Additional information:</b>                            | Not dangerous according to the above specifications. |
| · <b>UN "Model Regulation":</b>                                       | Void                                                 |

## SECTION 15: Regulatory information

### · 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

#### · Directive 2012/18/EU

· **Named dangerous substances - ANNEX I** None of the ingredients is listed.

#### · **DIRECTIVE 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment – Annex II**

None of the ingredients is listed.

#### · **REGULATION (EU) 2019/1148**

##### · **Annex I - RESTRICTED EXPLOSIVES PRECURSORS (Upper limit value for the purpose of licensing under Article 5(3))**

None of the ingredients is listed.

##### · **Annex II - REPORTABLE EXPLOSIVES PRECURSORS**

None of the ingredients is listed.

##### · **Regulation (EC) No 273/2004 on drug precursors**

None of the ingredients is listed.

##### · **Regulation (EC) No 111/2005 laying down rules for the monitoring of trade between the Community and third countries in drug precursors**

None of the ingredients is listed.

### · 15.2 Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

## SECTION 16: Other information

#### · **Date of previous version:** 14.02.2024

#### · **Abbreviations and acronyms:**

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative